Marvelous™ (MRVLS)

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References

- 1.Collignon PJ.Intravascular catheter associated sepsis: a common problem.
The Australian Study on Intravascular Catheter Associated Sepsis.Med J Aust 1994;161:374-8.
- 2.Heiselman D. Nosocomial bloodstream infections in critically ill patients.JAMA 1994;272:1819-20.
- 3.Norwood S, Ruby A, Civetta J, Cortes V. Catheter-related infections and associated septicemia.
Chest 1991;99:968-75.
- 4.Pittet P, Wendzel RP.Nosocomial bloodstream infections.Secular trends in mortality rates and contribution to total hospital deaths.Arch Intern Med 1995;155:1177-84.
- 5.O'Grady NP.Et al.Guidelines for the Prevention of Intravascular Catheter-Related Infections.CDC, 2002 Report.
- 6.Pearson ML.Intravascular device-related infections: an overview.Am J Infect Control 1996;24:262-93.
- 7.Ljungman P, Hagglund H, Bhjorkstrand B, Lonnqvist B, Ringden O. Perioperative teicoplanin for prevention of Gram-positive infections in neutropenic patients with indwelling central venous catheters: a randomized, controlled study.Support Care Cancer 1997;5:485-8.
- 8.Snydman DR, Pober BR, Murray SA, Gorbea HF, Majka JA, Perry LK.
Predictive value of surveillance skin cultures in total parenteral nutrition-related infection.
Prospective epidemiologic study using semi-quantitative cultures.Lancet 1982;1385-8.
- 9.Salzman MB, Isenberg HD, Shapiro JF, Lipsitz PJ, Rubin LG.A prospective study of the catheter hub as the portal of entry for microorganisms causing catheter-related sepsis in neonates.J Infect Dis 1993;167:487-90.
- 10.Peters G, Locci R, Pulverer G. Adherence and growth of coagulase-negative staphylococci in surfaces of intravenous catheters.J Infect Dis 1982;146:479-82.
- 11.Sitges-Serra A, Linares J, Perez JL, Jaurieta E, Lorente L. A randomized trial on the effect of tubing changes on hub contamination and catheter sepsis during parenteral nutrition.J Parenter Enteral Nutr 1985;9:322-5.
- 12.Linares J, Sitges-Serra A, Garau J, Perez JL, Martin R. Pathogenesis of catheter sepsis: a prospective study with quantitative and semiquantitative cultures of catheter hub and segments.J Clin Microbiol 1985;9:322-5.
- 13.deCicco M, Chiaradia V, Veronesi A, et al.Source and route of microbial colonization of parenteral nutrition catheters.Lancet 1989;2:1258-61.
- 14.Medical Anti-Bacterial Stopcock, Business plan, Elcam Medical, June 2004.
- 15.McArthur et al., Stopcock Contamination in an ICU, *The American Journal of Nursing*, Jan. 1975.
- 16.Blood residues analysis/ Igal Bar-Ilan Ph.D / Migal Ltd for Elcam Medical.
- 17.Bacterial growth characterization in saline solutions containing low concentration of blood in different time periods/ Aminolab Ltd.For Elcam Medical